


STR Authentication: Using the ATCC public STR Database

Brief Tutorial March 2016



Using the Public STR Database




STR Profile Background

- The ATCC STR database includes profile standards for all distributed cell lines.
- 8 loci are enough to authenticate a cell line for research and publication purposes.
- ATCC uses STR analysis to screen all human cell lines for authenticity and purity before distribution providing a true baseline for researchers using these cell lines.
- Comparing an STR profile using the ATCC database will provide a measureable relationship between the tested cells and accepted standard cell lines.
- Tumor and transformed cell lines are more prone to genetic drift which can accelerate with passage number, media content and other factors.
- STR profiles assume two alleles; the presence of more than two alleles in DNA from normal cells indicates genomic heterogeneity, which is typically equated with contamination or genetic instability. Some cell lines may have more than 2 alleles at a loci as they are generally not normal cells. See the example on the last page.

Cell authentication services are available from ATCC
www.atcc.org/str





Using the Public STR Database

SEARCH THE STR DATABASE

As part of our continuing efforts to characterize and authenticate the cell lines in the Cell Biology collection, ATCC has developed a comprehensive database of short tandem repeat (STR) DNA profiles for all of our human cell lines. [View our brief tutorial before starting.](#)

1. [STR Profiling Analysis](#)
2. [Matching Algorithm](#)
3. [Interrogating the Database](#)

Search by ATCC Number:

OR

Search by Amelogenin (AMEL) + at least 7 loci: Separate each allele entry with a comma (e.g., CSF1PO = 11, 12)

AMEL: ?

Plus

D5S818:

D13S317:

D7S820:

D16S539:

vWA:

TH01:

TPOX:

CSF1PO:

Go to the STR service landing page
or directly to the database

https://www.atcc.org/STR_Database.aspx?slp=1

There is a very simple registration
required to access the database.





Using the Public STR Database

There are 2 Choices

1. Select either search by ATCC number to obtain an STR profile standard

OR

2. Input a profile to run it against the ATCC public database

Search by ATCC Number:

OR

Search by Amelogenin (AMEL) + at least 7 loci: Separate each allele entry with a comma (e.g., CSF1PO = 11, 12)

AMEL: ?

Plus

D5S818:

D13S317:

D7S820:

D16S539:

vWA:

TH01:

TPOX:

CSF1PO:





1. Search by ATCC Number to Obtain an STR Profile Standard

Example:
ATCC catalog CCL-2 are
HeLa cells

Hit the search button

Search by ATCC Number:

OR

**Search by Amelogenin (AMEL) + at least 7 loci:
Separate each allele entry with a comma (e.g., CSF1PO = 11, 12)**

AMEL: ?

Plus

D5S818:

D13S317:

D7S820:

D16S539:

VWA:

TH01:

TPOX:

CSF1PO:





1. Search by ATCC Number to Obtain an STR Profile Standard

Example:

Find your CCL-2 HeLa cell
STR profile at the top of the
page

Export the data to
Excel

Showing 1 - 1 Of 1 PageSize: 100 ▼

Add to Cart	%Match	ATCC® Number	Designation	D5S818	D13S317	D7S820	D16S539	vWA	TH01	AMEL	TPOX	CSF1PO
<input type="checkbox"/>	100	CCL-2	HeLa	11,12	12,13,3	8,12	9,10	16,18	7	X	8,12	9,10

Disclaimer: Reference to this database and the data contained therein may be cited in publications, and ATCC encourages such citation or reference. While every reasonable effort has been made to assure the accuracy of these data, no warranty, express or implied, is made by ATCC as to their accuracy.

While ATCC has largely used the Promega PowerPlex® 1.2 System in the creation of these data and recommends that researchers wishing to produce data for comparison also use a Promega PowerPlex® System ATCC does not provide a general endorsement of this product or provide any warranty or representation regarding its quality or performance in the scientific community for the identification of human cell lines.

As in the past, when we find a misidentified cell line among our holdings (i.e., the DNA profile is similar or identical to that of an unrelated cell line), we will post a note on the Misidentified Cell Lines page of our website.

Search by ATCC Number:

OR

Search by Amelogenin (AMEL) + at least 7 loci:





2. Search by STR Profile

to Match Against Others in the Database

Input a profile to run against the ATCC public database

Search by Amelogenin (AMEL) + at least 7 loci:
Separate each allele entry with a comma (e.g., CSF1PO = 11, 12)

AMEL:



Plus

D5S818:

D13S317:

D7S820:

D16S539:

vWA:

TH01:

TPOX:

CSF1PO:

Matches >= 80%

Matches >= 56%

Clear





2. Search by STR Profile

to Match Against Others in the Database

Profile from a recent cancer research paper.

Separate alleles using a comma

For homozygous use a single number and NOT "10,10" for example.

Search by Amelogenin (AMEL) + at least 7 loci:
Separate each allele entry with a comma (e.g., CSF1PO = 11, 12)

AMEL: ?

Plus

D5S818:

D13S317:

D7S820:

D16S539:

vWA:

TH01:

TPOX:


CSF1PO:

To limit the number of results:


select either 80% match and higher
OR

56% match or higher





2. Search by STR Profile to Match Against Others in the Database



Results can take 30 seconds +/-

In this case there were no matches
below 80% for this profile

Showing 1 - 1 Of 1 Page Size: 100

Add to Cart	%Match	ATCC® Number	Designation	D5S818	D13S317	D7S820	D16S539	vWA	TH01	AMEL	TPOX	CSF1PO
<input type="checkbox"/>	80	CCL-171	MRC-5	11,12	11,14	10,11	9,11	15	8	X,Y	8	11,12

Disclaimer: Reference to this database and the data contained therein may be cited in publications, and ATCC encourages such citation or reference. While every reasonable effort has been made to assure the accuracy of these data, no warranty, express or implied, is made by ATCC as to their accuracy.

While ATCC has largely used the Promega PowerPlex® 1.2 System in the creation of these data and recommends that researchers wishing to produce data for comparison also use a Promega PowerPlex® System ATCC does not provide a general endorsement of this product or provide any warranty or representation regarding its quality or performance in the scientific community for the identification of human cell lines.

As in the past, when we find a misidentified cell line among our holdings (i.e., the DNA profile is similar or identical to that of an unrelated cell line), we will post a note on the Misidentified Cell Lines page of our website.

Search by ATCC Number:

OR

Search by Amelogenin (AMEL) + at least 7 loci: Separate each allele entry with a comma (e.g., CSF1PO = 11, 12)

AMEL: ?

Plus

D5S818:

D13S317:

D7S820:

D16S539:

vWA:

TH01:

TPOX:

CSF1PO:



2. Search by STR Profile Practice Profile

In this example you will find a
number of results only by selecting
 $\geq 80\%$

Search by ATCC Number:

OR

Search by Amelogenin (AMEL) + at least 7 loci: Separate each allele entry with a comma (e.g., CSF1PO = 11, 12)

AMEL:


Plus

D5S818:	<input type="text" value="12,13"/>
D13S317:	<input type="text" value="11,13,14"/>
D7S820:	<input type="text" value="11,12"/>
D16S539:	<input type="text" value="12,13"/>
vWA:	<input type="text" value="15,18"/>
TH01:	<input type="text" value="6"/>
TPOX:	<input type="text" value="9,11"/>
CSF1PO:	<input type="text" value="11"/>





2. Search by STR Profile to Match Against Others in the Database



Results can take 30 seconds +/-

In this case there were no matches
below 80% for this profile

Showing 1 - 3 Of 3 PageSize: 100

Add to Cart	%Match	ATCC# Number	Designation	D5S818	D13S317	D7S820	D16S539	VWA	TH01	AMEL	TPOX	CSF1PO
<input type="checkbox"/>	93.75	HTB-37	Caco-2	12,13	11,13,14	11,12	12,13	16,18	6	X	9,11	11
<input type="checkbox"/>	93.75	CRL-2102	C2BBE1	12,13	11,13,14	11,12	12,13	16,18	6	X	9,11	11
<input type="checkbox"/>	80	CRL-2868	HCC927	12	9	11,12	12	18	6	X	8	11

Disclaimer: Reference to this database and the data contained therein may be cited in publications, and ATCC encourages such citation or reference. While every reasonable effort has been made to assure the accuracy of these data, no warranty, express or implied, is made by ATCC as to their accuracy.

While ATCC has largely used the Promega PowerPlex® 1.2 System in the creation of these data and recommends that researchers wishing to produce data for comparison also use a Promega PowerPlex® System ATCC does not provide a general endorsement of this product or provide any warranty or representation regarding its quality or performance in the scientific community for the identification of human cell lines.

As in the past, when we find a misidentified cell line among our holdings (i.e., the DNA profile is similar or identical to that of an unrelated cell line), we will post a note on the Misidentified Cell Lines page of our website.

Search by ATCC Number:

OR

Search by Amelogenin (AMEL) + at least 7 loci: Separate each allele entry with a comma (e.g., CSF1PO = 11, 12)

AMEL:

Plus

D5S818:

D13S317:

D7S820:

D16S539:

VWA:


TH01:

TPOX:

CSF1PO:



2. Search by STR Profile to Match Against Others in the Database



Results can take 30 seconds +/-

In this case there were no matches 80% and above, so only the 56% and lower selection provided results.

Search by Amelogenin (AMEL) + at least 7 loci. Separate each allele entry with a comma (e.g., CSF1PO = 11, 12)

AMEL:

Plus:

D5S818:

D13S317:

D7S820:

D16S539:

vWA:

TH01:

TPOX:

CSF1PO:

Showing 1 - 100 Of 1139 Page Size:

Add to Cart	% Match	ATCC® Number	Designation	D5S818	D13S317	D7S820	D16S539	vWA	TH01	AMEL	TPOX	CSF1PO
<input type="checkbox"/>	31.25	CCL-2	HeLa	11,12	12,13.3	8,12	9,10	16,18	7	X	8,12	9,10
<input type="checkbox"/>	31.25	CCL-5	L-132	11,12	12,13.3	8,12	9,10	16,18	7	X	8,12	9,10
<input type="checkbox"/>	31.25	CCL-6	Intestine 407	11,12	12,13.3	8,12	9,10	16,18	7	X	8,12	9,10
<input type="checkbox"/>	35.71	CCL-13	Chang Liver	12	12,13.3	8,12	9,10	16,18	7	X	8,12	10
<input type="checkbox"/>	31.25	CCL-17	KB	11,12	12,13.3	8,12	9,10	16,18	7	X	8,12	9,10
<input type="checkbox"/>	35.71	CCL-21	AV3	11,12	13.3	12	9,10	16,18	7	X	8,12	9,10
<input type="checkbox"/>	31.25	CCL-23	HEp-2	11,12	12,13.3	8,12	9,10	16,18	7	X	8,12	9,10
<input type="checkbox"/>	31.25	CCL-23	HEp-2	11,12	12,13.3	8,12	9,10	16,18	7	X	8,12	9,10
<input type="checkbox"/>	33.33	CCL-25	WISH	11,12	13.3	8,12	9,10	16,18	7	X	8,12	9,10
<input type="checkbox"/>	52.94	CCL-30	RPMI 2650	12,13	11,12	8,11	11,12	16,18	6,8	X,Y	8	9,11
<input type="checkbox"/>	31.25	CCL-62	FL	11,12	12,13.3	8,12	9,10	16,18	7	X	8,12	9,10
<input type="checkbox"/>	42.86	CCL-75	WI 38	10	11	9,11	11,12	19,20	8,9.3	X	8	10,12
<input type="checkbox"/>	46.15	CCL-75.1	WI-38 VA13	10	11	9,11	11,12	19,20	9.3	X	8	10,12



2. Search by STR Profile to Match Against Others in the Database

<input type="checkbox"/>	37.5	CRL-1552	MOLT-3	12	12,13	8,10	11,13,14	17	6,8	X,Y	8	11,12
<input type="checkbox"/>	37.5	CRL-1740	LNcaP.FGC	11,12	10,12	9.1,10.3	11	16,18	9	X,Y	8,9	10,11
<input type="checkbox"/>	37.5	CRL-1740	LNcaP.FGC	11,12	10,12	9.1,10.3	11	16,18	9	X,Y	8,9	10,11
<input type="checkbox"/>	37.5	CRL-1740	LNcaP.FGC	11,12	10,12	9.1,10.3	11	16,18	9	X,Y	8,9	10,11
<input type="checkbox"/>	28.57	CRL-1598	A673	11,12	8,13	10,12	11	15,18	9.3	X	8	11,12
<input type="checkbox"/>	46.67	TIB-202	THP-1	11,12	13	10	11,12	16	8,9.3	X,Y	8,11	11,13

<< < 1 2 3 4 5 ... > >>

Sort by “% match” then check off the records you would like to keep.

Export your data to Excel

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O
1	% Match	Sample Count	Matches	Atcc Number	Designation	D55818	D13S317	D7S820	D16S539	vWA	TH01	AMEL	TPOX	CSF1PO	
2	31.25	16	5	CCL-2	HeLa	11,12	12,13.3	8,12	9,10	16,18	7	X	8,12	9,10	
3	31.25	16	5	CCL-5	L-132	11,12	12,13.3	8,12	9,10	16,18	7	X	8,12	9,10	
4	31.25	16	5	CCL-6	Intestine 407	11,12	12,13.3	8,12	9,10	16,18	7	X	8,12	9,10	
5	35.71	14	5	CCL-13	Chang Liver		12	12,13.3	8,12	9,10	16,18	7	X	8,12	10
6	31.25	16	5	CCL-17	KB	11,12	12,13.3	8,12	9,10	16,18	7	X	8,12	9,10	
7	35.71	14	5	CCL-21	AV3	11,12	13.3	12	9,10	16,18	7	X	8,12	9,10	
8	31.25	16	5	CCL-23	HEp-2	11,12	12,13.3	8,12	9,10	16,18	7	X	8,12	9,10	
9	31.25	16	5	CCL-23	HEp-2	11,12	12,13.3	8,12	9,10	16,18	7	X	8,12	9,10	
10	33.33	15	5	CCL-25	WISH	11,12	13.3	8,12	9,10	16,18	7	X	8,12	9,10	
11	52.94	17	9	CCL-30	RPMI 2650	12,13	11,12	8,11	11,12	16,18	6,8	X,Y	8	9,11	
12															
13															



Test your cells...
...Trust your data

Consider ATCC STR Cell Authentication