

**ATCC Medium: 2855 Tryptic Soy Agar/Broth with 100 mcg/ml Ampicillin**

**Agar Medium**

Tryptic Soy Agar (BD 236950).....40 g  
DI Water.....950 ml

Autoclave at 121°C and allow to cool down to 45-50°C. Aseptically add the following filter sterilized solution:

**Ampicillin Solution**

Ampicillin.....100 mg  
DI Water.....50 ml

Dispense as required. **Store at 4°C. Do not store longer than 7-10 days.**

**Broth Medium**

Tryptic Soy Broth (BD cat 211825).....30 g  
DI Water.....950 ml

Autoclave at 121°C and allow to cool down to 45-50°C. Aseptically add the following filter sterilized solution:

**Ampicillin Solution**

Ampicillin.....100 mg  
DI Water.....50 ml

Dispense as required. **Store at 4°C. Do not store longer than 7-10 days.**