



Product Sheet

Fusarium keratoplasticum (ATCC® 36031™)

Please read this FIRST



Storage Temp.
Frozen: -80°C or colder
Freeze-Dried: 2°C to 8°C
Live Culture: See Propagation Section



Biosafety Level
2

Intended Use

This product is intended for research use only. It is not intended for any animal or human therapeutic or diagnostic use.

Citation of Strain

If use of this culture results in a scientific publication, it should be cited in that manuscript in the following manner: *Fusarium keratoplasticum* (ATCC® 36031™)

American Type Culture Collection
PO Box 1549
Manassas, VA 20108 USA
www.atcc.org

800.638.6597 or 703.365.2700
Fax: 703.365.2750
Email: Tech@atcc.org

Or contact your local distributor

Description

Strain Designation: FIV/74 [NRRL 22641]

Deposited Name: *Fusarium solani* (Martius) Saccardo

Product Description: An ampoule containing viable cells (may include spores and mycelia) suspended in cryoprotectant.

Propagation

The information recommended in this section is to assist users in obtaining living culture(s) for their studies. The recommendation does not imply that the conditions or procedures provided below are optimum. Experienced researchers may initiate the growth of a culture in their own way.

ATCC® Medium 28: Emmons' modification of Sabouraud's agar

ATCC® Medium 200: YM agar or YM broth

ATCC® Medium 336: Potato dextrose agar (PDA)

Growth Conditions

Temperature: 24°C to 26°C

Atmosphere: Typical aerobic

Recommended Procedure

Frozen ampoules packed in dry ice should either be thawed immediately or stored in liquid nitrogen. If liquid nitrogen storage facilities are not available, frozen ampoules may be stored at or below -70°C for approximately one week. **Do not under any circumstance store frozen ampoules at refrigerator freezer temperatures (generally -20°C).** Storage of frozen material at this temperature will result in the death of the culture.

1. To thaw a frozen ampoule, place in a **25°C to 30°C** water bath, until just thawed (**approximately 5 minutes**). Immerse the ampoule just sufficient to cover the frozen material. Do not agitate the ampoule.
2. Immediately after thawing, wipe down ampoule with 70% ethanol and aseptically transfer at least 50 µL (or 2-3 agar cubes) of the content onto a plate or broth with medium recommended.
3. Incubate the inoculum/strain at the temperature and conditions recommended. Inspect for growth of the inoculum/strain regularly. The sign of viability is noticeable typically after 5-6 days of incubation. However, the time necessary for significant growth will vary from strain to strain.

Notes

This strain appeared to be a slow grower and less hardy than other *Fusarium* strains. Sporulation was also less and slower than other strains.

Additional, updated information on this product may be available on the ATCC® web site at www.atcc.org.

DNA Sequence

18S ribosomal RNA gene, partial sequence; internal transcribed spacer 1, 5.8S ribosomal RNA gene, and internal transcribed spacer 2, complete sequence; and 28S ribosomal RNA gene, partial sequence
GGTCTCGTTGGTGAACCGAGCGAGGGATCATTACCGAGTTATACAACCTCATCAACCTGTGAACATA
CCTAAAACGTTGCTTCGGCGGGAACAGACGGCCCCGTAACACGGGCGCCCGCCAGAGGACCCCC
TAACTCTGTTTCTATTATGTTTCTTCTGAGTAAAACAAGCAATAAATAAACTTCAACAACGGATCT
CTTGGCTCTGGCATCGATGAAGAACGCAGCGAAATGCGATAAGTAATGTGAATTGCAGAATTCAGTG
AATCATCGAATCTTTGAACGCACATTGCGCCCGCCAGTATTCTGGCGGGCATGCCTGTTTCGAGCGTCAT
TACAACCCTCAGGCCCGGGCCCTGGCGTTGGGGATCGGCGAGGCGCCCGCTGCGGGCACACGCCGTC
CCCCAAATACAGTGGCGGTCCCGCCCGAGCTTCCATTGCGTAGTAGCTAACACCTCGCAACTGGAGAG
CGGCGCGGCCACGCGTAAAACACCCAACTTCTGAATGTTGACCTCGAATCAGGTAGGAATACCCGCT
GAACTTAAGCATATCAATA

D1D2 region of the 28S ribosomal RNA gene

ATATCAATAAGCGGAGAAAAGAAACCAACAGGGATTGCCCCAGTAACGGCGAGTGAAGCGGCAAC
AGCTCAAATTTGAAATCTGGCTCTCGGGCCCGAGTTGTAATTTGTAGAGGATGCTTTTGGTGAGGTGCC
TTCCGAGTTCCTTGAACGGGACGCCATAGAGGGTGAGAGCCCGCTGTTGGTGGACACCGATCCTCTG
TAAAGTCCCTCGACGAGTCGAGTAGTTTGGGAATGCTGCTCTAAATGGGAGGTATATGCTTTCTAAAG
CTAAATACCGGCCAGAGACCGATAGCGCACAAAGTAGAGTGTATCGAAAAGATGAAAAGAACTTTGAAA
AGAGAGTTAAACAGTACGTGAAATGTTGAAAGGGAAGCGCTTGTGACCAGACTTGGGCTTGGTTGAT
CATCCGGGTTCTCCCGGTGCACTTCCGGCCAGGCCAGCATCAGTTCCGCTGGGGGACAAAGG

